

# REACH IN PLANT GROWTH

## PGC-40L2 Reach in Plant Growth Chamber



**Applications** This chamber product is frequently used for research applications such as lighting for vascular plants to facilitate standard plant production, plant pathology research and seedling germination and development.

Many other applications exist for this product. Please compare your own requirements to the specifications listed below.

**Controller** Percival's Intellus Ultra controller is capable of controlling temperature, humidity, CO<sub>2</sub> and lighting. The Intellus Ultra Control System is a single-board electronic solid-state design which includes a 10 key membrane keypad with LED indicators and a vacuum fluorescent display. Programs may be configured to run in real time or countdown (circadian) mode. Ramping and non-ramping program methods are available for each programming mode. Multiple programs can be linked to create complex environmental profiles. The Intellus Web Server (optional) allows for monitoring and controlling of the chamber via a web browser (requires Internet Explorer 6.0 +). This option allows for remote monitoring and programming of your chamber including alerts and current condition updates for up to five e-mail addresses. Please refer to [www.percival-scientific.com](http://www.percival-scientific.com) for additional information regarding the control system.

**Lighting System** There are two light fixtures per tier. The light fixtures yield up to 900  $\mu\text{moles}/\text{m}^2/\text{s}$  @ 6" from the lamps. The lamp provides a balanced spectrum for plant growth using T-5 fluorescent lamps plus extended life tungsten incandescent lamps. Programming and control of the lighting is done via Intellus real time controller. There are four levels of programming of fluorescent lighting and two levels of programming of incandescent lighting.

**Air Flow** Uniform air circulation moves across each shelf through air diffusers on the rear wall. Air travels along the entire back wall, over the shelves and returns to the ceiling fans through an opening between the light fixture and the door.

Temp Range (with all lights on)	Interior Space (volume)		Work Area		Maximum Growing Height		Exterior Dimensions in. (cm)			Light Intensity (6" from lamps unless otherwise noted)	# of Tiers
	° C	ft <sup>3</sup>	m <sup>3</sup>	ft <sup>2</sup>	m <sup>2</sup>	in.	cm	(W)	(D)		
10-44±1.0	110	3.1	40	3.72	30.5 *	142.2	100.5(255.3)	36(91.4) **	110.5(280.7)	900	2

\* 66.25" with center light fixtures and shelves removed.

\*\* Chamber depth does not include door thickness 1 5/8" and control panel which is 4".

# REACH IN PLANT GROWTH

## PGC-40L2 Reach In Plant Growth Chamber

**Cabinet Construction** All chambers are built in panel sections. Each section consists of 2" thick urethane insulation, metal interior and exterior surfaces, cam-type fasteners and vinyl gaskets. Interior and exterior are constructed of 22 gauge electro-zinc plated steel with the exception of the interior floor, which is stainless steel. The chamber cabinet is attached to an angle frame base which contains heavy duty swivel casters.

**Insulation** Woodless construction using foam-in-place 2" thick CFC free urethane insulation foam. This is an environmentally friendly foam with global warming potential (GWP) of 0.0 and ozone depletion potential (ODP) of 0.0.

**Doors** Two door openings each 29.25" x 57.25" (74.3 cm x 145.4 cm). A magnetic gasket provides tight seal to door frame.

**Interior Space** 110 ft<sup>3</sup> (3.1 m<sup>3</sup>) with a work area of 40 ft<sup>2</sup> (3.72 m<sup>2</sup>) provided on two tiers.

**Finish** Interior and exterior painted with highly reflective, environmentally friendly, high temperature baked white powder coating.

**Refrigeration** Self-contained water-cooled condensing unit with hot gas bypass system for continuous compressor operation, extended life and close temperature control. This continuous running condensing unit ensures precise temperature control by alternately cycling refrigerant and hot gas to the coil; this also prolongs the life of the compressor, and eliminates the risk of ice build up in the coil. Solenoid valves have an extended stem for quiet and long life operation. Heat rejection to the ambient by standard refrigeration system with a water-cooled self-contained condensing unit is 2500 BTU/hr. Optional outdoor all weather air-cooled condensing unit or self contained water-cooled condensing unit available upon request.

**Temperature Range** 4° - 44° C with all lights off and 10° - 44° C with all lights on (full fresh air) ± 1.0° C within the work area on a horizontal plane with all lights on.

**Temperature Safety Limit Controls** (Experiment Protection) Adjustable high and low temperature controls, audible alarms, and visual indicators are provided.

The controls shut down all the power to the chamber, and activates alarms. When the temperature returns to the normal range the system will automatically reset.

**Humidity Control (Optional)** Additive control of humidity in %RH through use of spray nozzles will maintain humidity levels up to 75% +/- 5% RH lights ON up to 95% +/- 5% lights OFF, between 15° and 30° C. Humidifier requires distilled or de-mineralized water. Optional dehumidification via independent coiling coil and reheat heaters will maintain humidity levels down to 40% RH between 15°C and 30° C.

**Options (most popular)** Advanced Intellus Control System (C9), Communications Software (C9+), Advanced Intellus with Touchscreen and Internet capabilities (C10), Spray nozzle humidifier with advanced RH sensor and some dehumidification via reheat heaters (H9), Dehumidification via independent cooling coil with reheat heaters and spray nozzle humidifier (H8), Ultrasonic Humidifier with advanced RH Sensor (H11), Dehumidification via dehumidifying coil with reheat heaters and Ultrasonic Humidifier (H12), Ultrasonic Humidifier with Electronic RH sensor (H14), CO<sub>2</sub> enrichment package, self contained air-cooled condensing unit, outdoor all weather air-cooled condensing unit or self-contained water-cooled condensing unit, Dry alarm contacts (S2), Closed loop dimmable lighting (Q22), Open loop dimmable lighting (Q23). Extended temperature ranges available. See other catalog sheets or consult factory for additional accessories.

**Convenience Receptacles** Two convenience receptacles provided inside chamber.

**Electrical Service Requirements** 120 - 208 VAC/3 phase/60 Hz, 4 wires plus ground - R.L.A. for standard chamber without any options is 24amps/leg (Provide MCA = 30 amps/leg minimum). Consult factory for electrical services when adding accessories to the chamber.



### ETA Associates

119 Foster Street, Bldg #6

Peabody, MA 01960

Tel: (978) 532-1330

Fax: (978) 532-7325

www.ETAassociates.com

eta@ETAassociates.com

